

DATA SHEET

Vet-PDA Occluder™

OCCLUSION PDA SYSTEM



PRODUCT: Implantable device to close a patent ductus arteriosus.

REFERENCES:

REF	Distal x Prox. Coil Diameter	Type	Coil Length (mm)	Implantation Catheter	Catheter Length
VETPDA000404	4 x 4 mm	FLEX	3.5	4 F	55 cm
VETPDA000504	5 x 4 mm	FLEX	3.5	4 F	55 cm
VETPDA000605	6 x 5 mm	FLEX	3.5	4 F	55 cm
VETPDA000706	7 x 6 mm	FLEX	4.5	4 F	55 cm

MATERIAL: Implantable Nitinol coil, medical grade stainless steel guide.

TECHNICAL CHARACTERISTICS: Vet-PDA Occluder™ has been designed for the safe and atraumatic closure of the PDA congenital heart defect “patent ductus arteriosus” of any type.

- Nitinol coil for ductus arteriosus up to 4 mm.
- Pre-assembled system.
- Distal reinforcement to prevent recoil.
- Release system with radiopaque markings.
- Flexible and adaptive design that allows precise and safe release of the coil.
- Repositionable and recoverable before release.
- Immediate closure.
- Single use only.
- Valid for all types of PDA.

PRESENTATION: Individual packaging

METHOD OF CONSERVATION: Normal temperature and humidity

STERILE: Yes

COUNTRY OF MANUFACTURE: Spain